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REFERENTIES

1. Battelino, T, Danne, T, Bergenstal, R M, Amiel, S A, Beck, R, et al. 2019. Clinical targets for continuous glucose monitoring data interpretation: Recommendations from the international consensus on time in range. *Diabetes Care*, 42, 1593-1603
2. DCCT/EDIC Study Research Group 2016. Intensive diabetes treatment and cardiovascular outcomes in type 1 diabetes: The DCCT/EDIC study 30-year follow-up. *Diabetes Care*, 39, 686-693.
3. Diabetes Control and Complications Trial Research Group 1993. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl Med J*, 329, 977-986.
4. Stratton, I M, Adler, A I, Neil H. A., Matthews D. R., Manley S. E., et al. 2000. Association of glycaemia with macrovascular and microvascular complications of type 2 diabetes (UKPDS 35): Prospective observational study. *BMJ*, 321, 405-412.
5. UK Prospective Diabetes Study (UKPDS) Group 1998. Intensive blood-glucose control with sulphonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet*, 352, 837-853.
6. Diabetes.Co.Uk. 2019. Diabetes and stress [Online]. Geraadpleegd op 23 maart 2021, van <https://www.diabetes.co.uk/diabetes-destress.html>
7. Vincze, G, Barner, J C & Lopez, D 2004. Factors associated with adherence to self-monitoring of blood glucose among persons with diabetes. *Diabetes Educ*, 30, 112-125.
8. Wagner, J, Malchoff, C & Abbott, G 2005. Invasiveness as a barrier to self-monitoring of blood glucose in diabetes. *Diabetes Technol Ther*, 7, 612-619.
9. Sagar, R C, Abbas, A & Ajjan, R 2019. Glucose monitoring in diabetes: From clinical studies to real-world practice. *Practical Diabetes*, 36, 57-62.
10. Hortensius, J, Kars, M C, Wierenga, W S, Kleefstra, N, Bilo, H J, et al. 2012. Perspectives of patients with type 1 or insulin-treated type 2 diabetes on self-monitoring of blood glucose: A qualitative study. *BMC Public Health*, 12, 167.
11. Davies, M J, D'alessio, D A, Fradkin, J, Kernan, W N, Mathieu, C, et al. 2018. Management of hyperglycaemia in type 2 diabetes, 2018. A consensus report by the American Diabetes Association (ADA) and the European Association for the Study of Diabetes (EASD). *Diabetologia* 2018;61(12):2461-2498.
12. Wright, L A & Hirsch, I B 2017. Metrics beyond hemoglobin A1c in diabetes management: Time in range, hypoglycemia, and other parameters. *Diabetes Technol Ther*, 19, S16-S26.
13. Nathan, D M, Kuenen, J, Borg, R, Zheng, H, Schoenfeld, D, et al. 2008. Translating the A1c assay into estimated average glucose values. *Diabetes Care*, 31, 1473-1478.
14. Beck, R W, Connor, C G, Mullen, D M, Wesley, D M & Bergenstal, R M 2017a. The fallacy of average: How using HbA(1c) alone to assess glycemic control can be misleading. *Diabetes Care*, 40, 994-999.
15. Rama Chandran, S, Tay, W L, Lye, W K, Lim, L L, Ratnasingam, J, et al. 2018. Beyond HbA1c: Comparing glycemic variability and glycemic indices in predicting hypoglycemia in type 1 and type 2 diabetes. *Diabetes Technol Ther*, 20, 353-362.
16. Aronson, R, Galystyan, G, Goldfracht, M, Sifri, S A, Elliot, L, et al. 2018. Direct and indirect health economic impact of hypoglycaemia in a global population of patients with insulin-treated diabetes. *Diabetes Research and Clinical Practice*, 138, 35-43.
17. Rossi, M C, Nicolucci, A, Ozzello, A, Gentile, S, Agliarolo, A, et al. 2019. Impact of severe and symptomatic hypoglycemia on quality of life and fear of hypoglycemia in type 1 and type 2 diabetes. Results of the Hypos-1 observational study. *Nutr Metab Cardiovasc Dis*, 29, 736-743
18. American Diabetes Association 2019. 7. Diabetes technology: Standards of medical care in diabetes-2019. *Diabetes Care*, 42, S71-S80.
19. Danne, T, Nimri, R, Battelino, T, Bergenstal, R M, Close, K L, et al. 2017. International consensus on use of continuous glucose monitoring. *Diabetes Care*, 40, 1631-1640.
20. Beck, R W, Riddlesworth, T, Ruedy, K, Ahmann, A, Bergenstal, R, et al. 2017b. Effect of continuous glucose monitoring on glycemic control in adults with type 1 diabetes using insulin injections: the DIAMOND randomized clinical trial. *JAMA*, 317, 371-378.
21. Beck, R W, Riddlesworth, T D, Ruedy, K, Ahmann, A, Haller, S, et al. 2017c. Continuous glucose monitoring versus usual care in patients with type 2 diabetes receiving multiple daily insulin injections: a randomized trial. *Ann Intern Med*, 167, 365-374.
22. Lind, M, Polonsky, W, Hirsch, I B, Heise, T, Bolinder, J, et al. 2017. Continuous glucose monitoring vs conventional therapy for glycemic control in adults with type 1 diabetes treated with multiple daily insulin injections: The GOLD randomized clinical trial. *JAMA*, 317, 379-387.
23. Van Beers, C A, Devries, J H, Kleijer, S J, Smits, M M, Geelhoed-Duijvestijn, P H, et al. 2016. Continuous glucose monitoring for patients with type 1 diabetes and impaired awareness of hypoglycaemia (IN CONTROL): a randomised, open-label, crossover trial. *Lancet Diabetes Endocrinol*, 4, 893-902.
24. Polonsky, W H, Hessler, D, Ruedy, K J & Beck, R W 2017. The impact of continuous glucose monitoring on markers of quality of life in adults with type 1 diabetes: Further findings from the DIAMOND randomized clinical trial. *Diabetes Care*, 40, 736-741.
25. Advani, A 2020. Positioning time in range in diabetes management. *Diabetologia*, 63, 242-252.
26. Vigersky, R A & McMahon, C 2019. The relationship of hemoglobin A1c to time-in-range in patients with diabetes. *Diabetes Technol Ther*, 21, 81-85.
27. Beck, R W, Bergenstal, R M, Riddlesworth, T D, Kollman, C, Li, Z, et al. 2019. Validation of time in range as an outcome measure for diabetes clinical trials. *Diabetes Care*, 42, 400-405.
28. Lu, J, Ma, X, Zhou, J, Zhang, L, Mo, Y, et al. 2018. Association of time in range, as assessed by continuous glucose monitoring, with diabetic retinopathy in type 2 diabetes. *Diabetes Care*, 41, 2370-2376.
29. Shrivastava, S R, Shrivastava, P S & Ramasamy, J 2013. Role of self-care in management of diabetes mellitus. *J Diabetes Metab Disord*, 12, 14.
30. Anderson, R J, Freedland, K E, Clouse, R E & Lustman, P J 2001. The prevalence of comorbid depression in adults with diabetes: A meta-analysis. *Diabetes Care*, 24, 1069-1078.
31. Gonzalez, J S, Safren, S A, Cagliero, E, Wexler, D J, Delahanty, L, et al. 2007. Depression, self-care, and medication adherence in type 2 diabetes: Relationships across the full range of symptom severity. *Diabetes Care*, 30, 2222-2227.
32. Runge, A S, Kennedy, L, Brown, A S, Dove, A E, Levine, B J, et al. 2018. Does time-in-range matter? Perspectives from people with diabetes on the success of current therapies and the drivers of improved outcomes. *Clin Diabetes*, 36, 112-119.
33. Diabetes.co.uk. 2019. Screening for Diabetes Complications [Online] Geraadpleegd op 30 juni 2021, van <https://www.diabetes.co.uk/diabetes-complications/screening-for-diabetic-complications.html>.